

Sotisu Medical Co., Ltd.

3-1-10-5F, Minami-Hakata Station, Hakata-ku, Fukuoka City, Fukuoka Prefecture

https://mytrex.jp https://usa.mytrex.jp

Designed in Japan. Made in China. ©2024 SOTISU MEDICAL CO., LTD.

Please read the manual carefully and use it correctly and safely. The removal of the battery cap is necessary. Please keep it properly after reading. For charge, warranty certificate attached.

Safety precautions. Please use it according to the prescribed methods, and violations are strictly prohibited. DANGEROUS. Prohibition: Do not use the following medical electronic equipment, please do not use...

MYTREX Customer Service. Business hours 9:00-17:30. [Japanese] Phone: 092-260-7366 Web: https://mytrex.jp/support/ [USA & Global] Email: Support.usa@mytrex.jp Web: https://usa.mytrex.jp

Warranty certificate. If there are any malfunctions after purchase, please contact our customer service. Warranty certificate Model: MT-EXP22B. The warranty period is 1 year under normal use.

WARNING. Do not use outside the designated areas and the following parts: Scalp surgery parts, Contactless body parts, Do not touch metal and other electrical products during use...

Quality Assurance. The warranty card is used according to the manual requirements under normal circumstances. This warranty certificate does not apply to accidents that occur during the warranty period due to the following circumstances.

WARNING. Prohibition: Do not use in the following conditions: In places with high humidity, In places with high temperature, In places with strong electromagnetic waves...

WARNING. Prohibition: Do not use together with a hair dryer or dry the main body and accessories with a hairdryer or other electric device. Do not use on people with skin conditions.

Product specifications. Power: 5V 1.1A, Battery: 1900mAh 3.7V, Main machine dimensions: 129x81x40mm, Material: Main machine: ABS, PC, Silicone, Attachments: Silicon, Carbon, Charger: ABS, TPU, USB cable: 1.2m, AC adapter: 1 unit.

Storage method. When you use, please make sure to wipe the water of the main unit and the attachments thoroughly. When not in use for a long time, please clean the main machine before storing.

WARNING. Prohibition: Do not use the AC adapter or USB cable if they are damaged or if the plug does not fit securely into the socket. Do not use in the bathroom.

WARNING. Do not use the machine and its attachments if they have abnormal odor or heat. Do not use the machine in a humid environment.

The EMS stimulation is strong and painful. Water remains inside after use. Do not use the AC adapter or USB cable if they are damaged.

Stop halfway. If part of the body or arms is trapped for more than 3 seconds during use, it will stop for safety reasons. The EMS stimulation is weak or I can't feel it.

WARNING. In case of any abnormal situation, please turn off the power supply of this machine. Before charging, please confirm that the charging terminals of the machine and the charger are free of moisture.

WARNING. For individuals undergoing dental treatment, to avoid receiving the treatment area when using this device, please consult your dentist. Do not use in the shower or bathtub.

Frequently Asked Questions. Power supply cannot turn on. Unable to charge. If the device does not start and make a "beep" sound when you press the power button...

5. Please use a dry cloth or towel to clean the water from the main machine and the attachments, and then keep them in a well-ventilated place to dry. 6. When the main body and attachment are dry, install the attachment to the main body.

WARNING. Do not use harsh irritants such as any friction creams. Do not use on people with skin conditions. Do not use in the bathroom.

WARNING. Do not use the device in the shower or bathtub. Do not use the device in high-temperature and high-humidity environments. Do not use the device in a humid environment.

4. Use the attached cleaning brush to remove stains on all parts and rinse with water. Clean the movable parts and the attachments of the main machine as shown in Figure 2.3.

Cleaning methods. 1. Make sure to turn off the power first. 2. Please wash the hair and stains on the main machine with water. 3. Remove the attachments from the main machine.

NOTICE. Always clean the machine after using disinfectant, soap treatments, lotions and gels. Please be sure to conduct a test before initial use. Do not use in the shower or bathtub.

NOTICE. For accessories (necklaces, earrings and etc.), please remove them before use. Please be sure to remove make-up when using on the face. Handle with care.

8. Before care the scalp, please wash it again with shampoo (conditioner and etc.). This product (main unit) contains a lithium-ion battery. Please read the manual carefully.

Body Care. Please use the scalp attachment on the body, arms, waist, thighs, legs and other parts. Do not use on people with skin conditions. Do not use on non-designated parts.

Name and function of each part. Main machine: Power button, EMS intensity button, Blue LED light, Green LED light, Red LED light, Blue LED light. Battery level indicator.

Attaching the device. 1. Make sure to turn off the power first. 2. Please wash the hair and stains on the main machine with water. 3. Remove the attachments from the main machine.

Facial care. Please be sure to use the face attachment. Please move slowly in the direction of the arrow for face, head, chest and neck care. Prohibited parts: Nose, ears, lips, eyelids, hair, eyebrows, forehead, throat.

Before use (Charging method). For initial use and when the battery is low, please turn off the power of the machine before charging. Please use the scalpel attachment.

Attaching. Care for the front and side of the head. Care for the back of the head. Use the same order and method as Step 1, start from the top of the head to the side.

Methods of use. The recommended water temperature is 5-42 degrees Celsius. 1. Select the attachments to install on the main machine. 2. Please fully moisturize the scalp and skin.

4. After charging, unplug the AC adapter from the socket. About charging: Do not use for a long time if a significantly longer and longer after.

Read carefully for the first time using. About Safety Device: If any part of the body or other objects is trapped or caught for more than 3 seconds. About attachments: Please use the scalp and face attachments alternately for various purposes.

3. Please confirm that there is no pain or hot flashes on the skin within 24 hours. If the product may cause skin redness due to different levels of stimulation after use.

Skin pre-test. When using this product for the first time, in order to prevent skin problems, please conduct a pre-test first. Please use the scalpel attachment.

Attachment notes. In order to average the EMS (current stimulation), tone is used on the attachments. Please be aware that black stains may adhere to the skin and face.

About EMS. EMS only relies on electrical stimulation of muscles, it's important to start at EMS level 1 (LED3) right away when using on the neck.

On rollers, + electrode and - electrode, as shown in the right picture. Although this product can be used on the scalp, face and body, the sensation of electrical stimulation varies greatly.

EMS works through water or moisture conduction. EMS stimulation is difficult to feel when the skin is dry or has hair on it. It is important to start at EMS level 1 (LED3) right away when using on the neck.

If there is no problem with the redness becoming lighter or even if the symptoms do not recover, please check with a doctor. This product may cause skin redness due to different levels of stimulation after use.