

Sliding Zero Gravity

Massage Chair

MC-918

User Manual



CONTENTS

| | |
|---------------------------------------|--------------|
| <u>IMPORTANT SAFETY INSTRUCTIONS</u> | 02-04 |
| <u>SPECIFICATIONS</u> | 05 |
| <u>STRUCTURE DESCRIPTION</u> | 06-07 |
| <u>SETTING UP THE MASSAGE LOUNGER</u> | 08 |
| <u>CHAIR INSTALLATION</u> | 09-10 |
| <u>BEFORE SITTING DOWN</u> | 11-12 |
| <u>REMOTE CONTROL</u> | 13 |
| <u>LCD DISPLAY</u> | 14 |
| <u>OPERATING INSTRUCTIAN</u> | 15-16 |
| <u>AFTER COMPLETING THE MASSAGE</u> | 17-18 |
| <u>OPERATION METHOD</u> | 19 |
| <u>CLEANING AND MAINTENANCE</u> | 20 |
| <u>WARRANTY CLAIMS</u> | 21 |

WARRANTY CLAIMS

1. **Warranty is not transferable;**
2. **This warranty is effective only if the product is purchased from Mas-Agee;**
3. **Proof of purchase(original receipt)is required for all warranty repairs;**
4. **All implied Warranties, including but ont limited to these implied warranties of fitness and merchantability, are limited to 1 year from the original purchase date;**
5. **Damage due to misuse, improper treatment and nuauthorised modifications and repairs are not covered by this warranty;**
6. **Warranty does not cover accessories and attachments that do not belong to the massage chair;**
7. **Warranty is not applicable to rental, business, commercial, institutional, or other monresidential users;**
8. **All services covered by this warranty must be approved by Mas-Agee and repairs are to be carried out by authorised technicians only.**
9. **Mas-Agee reserves theright to make substitution in lieu of repair or replacement, if replacement parts of defective materials are not available.**

IMPORTANT SAFETY INSTRUCTIONS

Basic precautions should always be followed when using an electrical unit, including the following.

⚠ DANGER – To reduce the risk of electric shock:

1. Always unplug the unit immediately after use.
2. Always unplug the unit before cleaning.

⚠ WARNING – To reduce the risk of burns, fire, electric shock, or injury to persons:

1. The unit should never be left unattended when plugged in. Always unplug the unit when not in use.
2. Close and continuous supervision is necessary when the unit is used by people who are limited in their ability tomove or communicate.
3. Only use the unit for the purpose described in these instructions.
4. Do not stand on the unit.
5. Do not use any accessories other than those recommended by the manufacturer.
6. This unit should not be used by children. Keep children away from the unit. Remove key when not in use and store out of reach of children.
7. Never operate the unit if it has a damaged power cord or plug. Return the unit to the nearest authorized service center if it is not working properly, if it has been dropped or damaged, or if it has been immersed in water.
8. Always raise the comfort cushion to check that the fabric of the unit itself has not been ripped before using the unit. Please also check other areas to ensure that the fabric has not been ripped. No matter how small the rip, if you find one then immediately stop using the unit, unplug the power cord and have the unit repaired at an authorized service center.
9. Keep power cords away from heated surfaces. Do not carry the unit by the power cord. Take care that pins,garbage and moisture do not come into contact with the plug. (Wipe off the dust using a dry cloth.)
10. The unit should not be used on top of heating appliances, such as electric carpets, etc. Do not operate the unit under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
11. Never block the air openings during operation and make sure that all air openings are kept free of lint and hair,etc. Do not drop or insert anything into any of the openings of the unit.
12. The unit should not be used in bathrooms or other damp or humid places as this may cause electric shock or the unit to malfunction. Water, etc. should not be spilled on the controller.
13. Do not use the unit outdoors.
14. Do not use the unit in places where aerosol spray products are in use or where oxygen is being used.
15. When unplugging the unit, turn off all controls before removing the plug from power outlet. A transformer should not be used as this may cause an electric shock or the unit to malfunction.
16. Connect this unit to a properly grounded outlet only. See Grounding Instructions.
17. The unit contains no user-serviceable parts and opening or disassembling any parts of it may cause an electric shock or injury.

SAVE THESE INSTRUCTIONS

IMPORTANT SAFETY INSTRUCTIONS

⚠ WARNING

- Persons with any of the conditions listed below or who are undergoing medical treatment should consult with their physicians before using the unit:
 - Pregnant women, people who are ill, in poor physical shape.
 - People suffering from back, neck, shoulder, or hip pain.
 - People who have a pacemaker, defibrillator, or other personal medical devices prone to electronic interference.
 - People suffering from heart disease.
 - People who are prohibited from receiving massage by a physician due to thrombosis or aneurism, acute varicose veins, or other circulatory disorders.
 - People with irregular curvatures of the spine.
 - People who have a loss of heat sensation.
- Make sure that no part of your body is between the massage heads as their squeezing action may cause injury. Your fingers, feet and head should not be put in the gap between the backrest and the seat, the backrest and the armrest, the seat and the cover under the armrest, the legrest and the footrest, the back cover and the leg cover, and the back side of the legrest.
- The unit should not be used against bare skin. While thin clothing may increase effectiveness, exposing the skin directly to the unit may irritate the skin.
Avoid wearing anything hard on your head like a hair accessory, etc. when using this unit. Excessively strong massage action should be avoided on the back of the neck.
- Avoid massaging any one part of your body for more than 5 minutes at a time.
Excessive massaging can overstimulate the muscles and nerves and result in an adverse effect or injury. The unit should not be used for more than a total of 16 minutes each session to avoid excessive massaging.
- While using the unit, if you start feeling sick or if the massage seems painful, stop using it immediately.
- You should not fall asleep while using the unit or use the unit after drinking alcohol.
- Avoid sitting on the unit, or sitting with your legs resting on the headrest, when the seat back is in a reclined position. Sitting on the legrest may cause the unit to fall over.
- Be sure to start on the gentle massage program if using the unit for the first time.
- If the unit functions abnormally, immediately turn off the power and have the unit checked by an authorized service center.
- Avoid attempting to repair the unit yourself.
- Make sure there are no obstacles behind the unit before reclining it. If the back hits a wall or pillar, the unit may malfunction. Be sure to have ample space behind the unit.
- After each massage, slide the power switch to the "off" position and turn to the "OFF" position.
- Be sure to fix a pillow on the back cushion for safe use.
Failure to do so may result in entangling your hair in the fastener, harm to your skin, or damage to your clothes.
- Confirm that the shoulder position is correct when body scanning is performed. If it is not properly aligned, use the POSITION button to adjust the shoulder position, or stop the massage and perform the body scanning again.
(For Automatic courses or Manual Operation)

CLEANING AND MAINTENANCE

⚠ Caution

- After finish, remember to pluck the plug.
- Before clean, remember to pluck the plug.



Avoid cause fire & electric shock.

Areas covered with synthetic leather

Wipe these areas with a soft, dry cloth.

- When using the commercially available leather-care products (wipes), follow its instructions.
- If the synthetic leather is particularly dirty, wipe it in the following manner.
 - Soak a soft cloth in water or a 3% to 5% solution of mild detergent and wring it out thoroughly.
 - Dab the surface with the moistened cloth.
 - Rinse the cloth in water, wring it out thoroughly, and then wipe any remaining detergent from the surface.
 - Wipe the surfaces with a soft, dry cloth.
 - Allow to dry naturally.
- If it is difficult to remove dirt, soak a commercially available Melamine foam sponge in mild detergent and wipe the unit with it.
- Do not use a hair dryer to dry the surface more quickly.
- Since colors may transfer from clothing to the synthetic leather surface, be careful when using the unit with clothing such as jeans or colored clothes.
- Do not allow these areas to come in contact with plastic for extended periods of time as this may lead to discoloration.
- Never use chemicals such as thinner, benzene, alcohol, etc.

⚠ Caution

No petrol, diluter, spirit to wipe,
No insecticide




Avoid cause fire & electric shock.



Side plate

- Wipe the unit with a cloth moistened with mild detergent and thoroughly wrung out.
Never use chemicals such as thinner, benzene, alcohol, etc.
- Wipe the unit with a cloth moistened with plain water and thoroughly wrung out.
Make sure to wring tightly before wiping the controller.
- Allow the unit to dry naturally.



OPERATION METHOD

Using trundle to move massage chair

- Lift up the Calves Rest.
- Then recline the chair to approximately to 45 degree, so that you can easily move the chair around the house with a slight pull or push.

Figure 1





Figure 1



It maybe harm to the floor, therefore , it is recommended that put some cushion on the floor.

- * If no anchor-hold chair leading drop out, it may harm to the floor.
- * Do not catch armrests when moving, as it easy to broken.

⚠ Caution

Loading the chair on the floor

- Loading on the floor Slowly load, be careful your foot.
- Anchor-hold the calve, after it load on the floor.
If calve not stretch completely, it will automatic return to initial position.

IMPORTANT SAFETY INSTRUCTIONS

15. Sit on the unit after confirming there is no foreign matter between the parts of the unit. (Confirm that there are no foreign objects stuck in the back cushion, legrest or hand/arm massage section.) Sit down after checking that the legrest does not rotate. Failure to observe this precaution may result in accident or injury.
16. Remove hard objects, such as false fingernails, wristwatches or rings, when massaging hands or arms. Do not massage with any hard objects in your pants pocket.
17. Store the attachment screws out of reach of children. Accidental ingestion may occur.
18. Do not stand, jump on, or place objects on the seat until the air is completely extracted for the air massage.
19. Do not pull on the legrest rotation lock release lever while sitting in the chair with your feet on the legrest. Accidents or injuries may occur.
20. Do not drag or push the unit in an installed state. Do not move the unit using the wheels on the floor. The floor may be damaged.
21. Do not hold on the legrest leg/sole massage section when moving the unit. If this is not observed, the section may slide as you move the unit, and injuries may occur.

GROUNDING INSTRUCTIONS

The unit must be properly grounded. If it should malfunction or break down, grounding provides a path of least resistance for electric current in order to reduce the risk of electric shock. This unit is equipped with a cord with an equipment-grounding conductor and a grounding plug. The plug must be plugged into an appropriate outlet that has been properly installed in accordance with all local codes and ordinances.

⚠ DANGER

Improper connection of the equipment-grounding conductor may lead to electric shock. Check with a qualified electrician if you are in doubt as to whether or not the unit is properly grounded. Do not modify the plug provided with the unit if it will not fit the power outlet. Instead, have a proper outlet installed by a qualified electrician.

This unit has a grounding plug. Make sure that the product is connected to a power outlet with the same configuration as the plug. Do not use an adaptor with this unit.

SPECIFICATIONS

| | |
|-----------------------------|---|
| Product Name | Sliding Zero Gravity Massage Chair |
| Model number of the product | MC-918 |
| Power Supply | AC220V 50/60Hz |
| Power Consumption | 120W |
| Maximum operating period | 30 minutes |
| Up or Down massage | about 130cm |
| Left and right massage | The length of massage wheels (Adjustable) neck、shoulder、back、waist: around7cm-20cm when it working, the length of massage wheels (include the width of massage wheels) : about7cm-20cm |
| Airbags adjustment | five-gear |
| Massage intensity | five-gear |
| Dimensions(H × W × D) | Main chair carton:1590*775*1070(cm) |
| Weight in box | N.W.:95KGS G.W.:112KGS |
| Length of the wire | Controller wire:1.7m Power supply wire:2.2m |
| Usage condition | Environment temperature:10°-40° Contrasting humidity:30-85RH |
| Storage condition | Storage temperature:20°-60° Storage humidity level:30-85RH |
| Maximum user weight | Approx. 120 KG |
| Minimum distance from wall | 2 cm |
| Using benefits | Increasing blood circulation Relieving muscular fatigue |


For those who weighing 100 kg or more, the product operations may make more noise and the fabric cover may wear faster.

AFTER COMPLETING THE MASSAGE


Disconnecting the power

1. Slide the power switch to the "O" position.

- Massage balls return to initial position .
- Calve auto decrease, back auto increase. (After massage balls return, the back will move.)



2. Remove the power plug from the electrical outlet.



AFTER COMPLETING THE MASSAGE

Timer function


The massage will complete automatically approximately 15 minutes after the start of operations.

- Sole massage section will retract, and the legrest will be lowered automatically so you can stand up safely.

Sole massage section will not retract all the way when the legrest is extended out, since the sole massage section will hit the floor.

Raise your legs to retract all the way.


- Backrest will not rise automatically.
- The massage heads move to the retracted position. Backrest will rise by pressing  button twice.



Interrupting the massage

Return the sole massage section to the original position and then press the  button.

- The massage heads move to the retracted position.
- The legrest lowers and the backrest raises automatically. (The backrest will move after the massage heads retract.)



Lift your soles slightly and legrest returns.


If you experience any problems

Press the  button.

- All operations will stop immediately.
- Get off the unit carefully to prevent the unit tipping over.

Returning to normal chair position


Hang the comfort cushion on the backrest.



STRUCTURE DESCRIPTION

Main unit


Front



- | | |
|---|--|
|  Pillow Pad |  Legrest |
|  Bluetooth Speaker |  Footrest |
|  Backrest Pad |  Backrest |
|  Arm air pressure |  LED Logo |
|  Controller |  Armrest |
|  Seat Pad |  Base |

STRUCTURE DESCRIPTION

Massage position(Enjoy the airbags full-wrap massage)




OPERATING INSTRUCTION

9. Massage chair angle adjusting control key

- 
The leg length adjustment control button
 Through this two key to control the angle of backrest and legs , according to a " " button, massage Angle will be adjusted. When the adjustment to the appropriate location, press the button of " "once again and positioning.

- 
Leg Angle adjusting control keys
 Through this two keys to control the leg Angle, according to a " " button, the leg frame Angle will be adjusted. When the adjustment to the appropriate location, press the button of " "positioning.



Caution

Description: when pressing the legs to lie down, fall and sit after the control key agreement 1 seconds or so, the back and leg angle adjustment.

10. Zero gravity adjustment

According to the " " button, the backrest and leg Angle adjustment to the zero gravity position, press the " "button, massage chair back to the initial position;

11. The leg length adjustment control button

Through this two control key can adjust the length of the leg massage, press a " " button, the calf length will be adjusted. When adjust the appropriate location, press again " "button to positioning.

Note:

According to your own needs, please adjust the angle or strength of the massage chair or cushion to reach the best massaging effect.

OPERATING INSTRUCTION

Switch on , press the power switch which under massage seat to "I "(open) position, the chair in reset state , The back turn to highest point , leg pad to lowest , the chair enter in ready mode. Below is massage controller's Button instruction :

⚠ Caution

- The unit must not be used for more than a total of 15 minutes.
 - Avoid massaging any one part of your body for more than 5 minutes at a time.
- Failure to observe this may result in an adverse effect or injury.

1. Press open the power.

The Start / Stop Control Button

Press once on this button and the backlight of the remote control will start blinking: the chair is in idle mode. Press any massage function button to begin massage function on the chair. Press this button again to stop all massage functions, the massage rollers will return to its original position.

2. The Pause Button.

When the chair is in massage mode and the user wants to temporarily pause it, press the pause button and all massage functions will stop.

There are five auto programs as follows:

3. **Stretch** (Focus on the foot stretching massage)
4. **Neck** (Focus on the neck massage)
5. **Man** (Suit for people who like high intensity of massage)
6. **Woman** (Suit for people who like low intensity of massage)
7. **Waist** (Focus on the waist massage)
8. **Sleep** (After 15mins massage,the chair stop at the zero gravity angle,power off,the chair reset)


⚠ Caution

- During the massage state when choose another mode .it will change to another function.
- In the automatic procedure, foot massage will start automatically.
- You can also change the strength of the leg massage or make it stop working.
- Press the automatic program then does massage automatic programs.

SETTING UP THE MASSAGE LOUNGER

Confirm the place

- Make sure there is enough room for massage chair incline to work.
- Do not expose massage chair in direct sunlight or high temperature environments (such as in the front of heating equipment). That may cause color fading or hardening of PU leather.
- Pave the mat to avoid damage by the floor.
- When paving the mat, the size of mat can cover the contact range.
- Choose the safe power socket
- Adopting the design of one-key forward sliding, the massage chair can be completely opened at a distance of 2cm against the wall.In the process of flipping,forward sliding can be reached to 26cm,the zero gravity angle can be adjusted to 158 degrees.



Take out parts and install the chair


- The chair is heavy(about 80kgs).Move the chair carefully to avoid person injured or damage of chair.

Put the chair on the floor

- Lay the chair slowly and watch out your feet.
- The floor maybe damaged, so we advise you to pave the mat or isolation pad.


CHAIR INSTALLATION

Confirm that the main unit and the accessories are in the box.



LCD DISPLAY

- * According to the area of the machine bought different, the LCD control panel display language may be different.
- * Remote control screen by capacitive touch screen plus break code LCD screen.



- 1 Message remaining time display**
- 2 Pause button display**
 - According to suspend tag will shine.
- 3 Massage position display**
 - The currently selected gasbag massage on operation;
 - When lit means gasbag massage position. Figure 2
 - When lit means the parts of body massage. Figure 1




Figure 1 Figure 2

- 9 Massage technique display**
 - Massage process, display the current icon.
 - The automatic program control, massage technique icon lit.



- 10 Back strength display (in kneading massage mode)**
 - Massage strength points 5 levels display
- 11 Airbag strength display**
 - Airbag strength points 5 levels display

- 4 Automatic display program**
 - Open air bag massage from head to toe
- 5 Back massage range**
- 6 Heating massage display**
- 7 Back fix point massage display**
- 8 Bluetooth display**
 - The Bluetooth Identification name is XY_BT.


REMOTE CONTROL

*Display language on the controller may differ depending on the region the unit was purchased.




- 1 On/Off button**
 - Switches the controller on and off.
- 2 Stop button**
 - Press this button to stop the function.
- 3 Stretch**
 - Open or close the “Stretch” program.
- 4 Neck**
 - Open or close the “Neck” program.
- 5 Man**
 - Open or close the “Man” program.
- 6 Woman**
 - Open or close the “Woman” program.
- 7 Waist**
 - Open or close the “Waist” program.
- 8 Sleep**
 - Open or close the “Sleep” program.
- 9 Angle Adjustment button**
 - Calves/Foot/Back angle adjustment button
- 10 “Zero gravity” button**
- 11 Foot flexible adjustment button**
- 12 Touch screen display**

CHAIR INSTALLATION





1. In the right side of the armrest front find remote control bracket mounting hole, lock screw of the specifications of the screw is:M6x45 2pcs;
2. The cushion of the trachea and vibration line connected in color, the zipper pull well;
3. Turn the legrest over from the seat and connect with the seat connecting rod, insert the pin.
4. Zip up the cushion.



CONNECT THE POWER CORD TO THE MACHINE

- Check power input socket in power box and whether power line interface is dirty or damaged;
- Put the power plug into the power input socket;
- Make sure the plug insert fully and safely;

Put power plug into a power outlet




BEFORE SITTING DOWN

Check the following before using the unit

1. Check the surrounding area.

Make sure there are no objects, people, or pets in the vicinity of the unit.




⚠ Caution


Avoid cause machine fault, accident or injury:

- Please don't put your hand or foot in the base before adjusting the height of footrest.
- Make sure there is no obstacle between footrest and also in front of footrest and underneath the footrest.
- Make sure there is no obstacle in the back of the backrest when adjusting the angle of backrest.
- Please don't sit or stand on the elevated footrest or backrest.

2. Check the power cord and the power plug



Be sure to push the connector all the way into the appliance inlet.




⚠ Caution

- The product must be grounded. If the machine fails dysfunction, the grounding will provide for the current path of least resistance, to reduce the risk of electric shock.
- If the equipment power plug is not properly connected, it will cause electrical hazard. If you are in doubt of how to correct the product, see a professional for guidance.

BEFORE SITTING DOWN

3. Switching on the power supply

1. Insert the power plug into an electrical outlet.
2. Slide the power switch to the "I" position.




4. Check the unit

1. Make sure there are no tears or rips in the fabric.

Shoulder/Arm massage position

- Make sure there are no tears or rips in the fabric.
- Make sure there is no foreign matter between the unit.
- Make sure the legrest is in lowest position.



⚠ Caution


Note: If you sit when the legrest is in lift position, it may make the hip to meet the connection between the legrest and the seat, or it will be hurt due to the chair turnover.

2. Make sure there is no foreign matter between the unit.
3. Make sure the massage heads are in the retracted position.

When the massage heads are not in the retracted position

Press the  button twice to return the massage heads to the retracted position.

Retracted position



4. Make sure the legrest is in the fully lowered position.

When the legrest is not fully lowered

The legrest will lower by pressing on  twice.

Check